RetinaStation

User Manual Manual del usuario Manuel de l'utilisateur Manuale utente Benutzerhandbuch

CE
0120

Users are required to read through the accompanying documents before using this instrument.

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1 Introduction

1.1 Outline

The Nikon non-mydratic fundus camera RetinaStation is a fast and easy to use auto 3D tracking retinal imaging system. RetinaStation is designed to provide images of the eye as an aid to clinicians in the diagnosis of diabetic retinopathy, AMD, glaucoma and other retinal diseases.

1.2 Intended Use

RetinaStation provides non-mydratic color retina and external images of the eye as an aid to clinicians in the evaluation and diagnosis of eye disease. RetinaStation, which uses NIR LED as illumination when aligning to the retina of a patient's eyes, only provides images and does not provide any diagnostic or pathological analysis or conclusions in any way. Furthermore, RetinaStation is not intended to be used as the sole diagnostic aid in disease identification, classification or management.

1.3 Proper Use of the Instrument

1. RetinaStation is a medical instrument; it must be operated by properly trained and qualified person(s) only. Its operation should be supervised by a physician.
2. Please be sure to read the user manual and understand the safety precautions before operating this instrument.
3. Always enter patient information first.
4. Prepare the patient contact surfaces (forehead rest and chin rest) according to the cleaning method described in this manual.
5. Immediately turn off the power switch of this instrument and disconnect the power cable if unusual problems arise.
6. Clean the ocular lens frequently to ensure good image quality.
7. Properly adjust the height of the motorized adjustable table to ensure the comfort of the patient during examination.
8. Align the patient's eye position to the canthus indicator mark on the chin and forehead rest assembly.
9. Dim the room lights to allow natural dilation of the patient's pupil and provide comfortable viewing of the fixation target without glare.
10. System functionality should be inspected once a year or whenever any repairs are made.

2 Safety Information

2.1 Displays for Using the Instrument Safely

<u>Display</u>	<u>Meaning</u>
 WARNING	A WARNING indicates the presence of a hazard that could result in severe personal injury.
 CAUTION	A CAUTION indicates the presence of a hazard that could result in minor injury.
NOTE	A NOTE provides useful information for operation that is important.

 WARNING	Accessory equipment that connects to the digital interfaces must be certified according to the respective IEC standards (e.g., IEC 60950 for laptops or IEC 60601-1 for medical equipment). Furthermore, all configurations shall comply with the system standards IEC 60601-1-1 and IEC 60601-1:2005. Any person who connects/installs devices to/on the system is responsible for verifying compliance. If in doubt, consult your local Nikon representative or distributor.
 WARNING	To avoid the risk of electric shock, this equipment must only be connected to a supply mains with a protective earth.
 WARNING	Do not modify this instrument without the authorization of the manufacturer.
 WARNING	THE RETINASTATION CANNOT REPLACE CLINICAL JUDGEMENT AND IS INTENDED TO BE USED ONLY IN CONJUNCTION WITH OTHER CLINICAL TOOLS CONSIDERED TO BE THE STANDARD OF CARE FOR MEASUREMENT AND DIAGNOSIS OF THE EYE.
 WARNING	The RetinaStation is a medical instrument. The software and hardware has been designed in accordance with U.S., European and other international medical instrument design and manufacturing standards. Unauthorized modification of the RetinaStation software or hardware, or any addition or deletion of any application in any way can jeopardize the safety of operators and patients, the performance of the instrument, and the integrity of patient data. <u>Any changes, additions or deletions to factory installed applications, operating system or modifications to hardware in any manner VOIDS the warranty completely.</u>
 WARNING	RetinaStation is not intended for home use and may not be stored or operated in environment conditions other than those prescribed (see "Specifications").

 WARNING	<p>Phototoxicity Prolonged intense light exposure can damage the retina. Therefore, the use of the device for ocular examination should not be unnecessarily prolonged, and the brightness setting should not exceed what is needed to provide clear visualization of the target structures.</p> <p>The retinal exposure dose for a photochemical hazard is a product of the radiance and the exposure time. If the value of radiance is halved, twice the time is needed to reach the maximum exposure limit.</p>
 WARNING	<p>Do not obstruct the mains power switch or position the equipment where the connection to the mains line can be accidentally disconnected.</p>
 WARNING	<p>The instrument is not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide.</p>
 WARNING	<p>The RetinaStation has no special protection against the harmful ingress of water or other liquids (classified IPX0). To avoid damaging the instrument and causing a safety hazard, cleaning solutions, including water, should not be directly applied to the instrument. Using a dampened cloth (without dripping) is a good method for cleaning the exterior surface of the enclosure.</p>
 WARNING	<p>The patient cannot touch any electrical device that is not powered by RetinaStation with any part of his or her body while being examined. In addition, the RetinaStation operator must not attempt to touch the patient and any electrical device that is not powered by RetinaStation while examining the patient. Failure to do so could result in electrical shock to the patient and/or operator.</p>
 WARNING	<p>Do not connect any unspecified devices to the instrument. Doing so may result in fire or electric shock. For details regarding purchasing accessories, please contact your local Nikon representative or distributor.</p>


 CAUTION	<p>Be sure to hold the bottom of the base when moving RetinaStation.</p>
 CAUTION	<p>Do not operate the LCD monitor with wet hands or hard objects. Doing so may damage the LCD monitor.</p>
 CAUTION	<p>When adjusting the chin rest, be careful not to pinch the patient's hands.</p>
 CAUTION	<p>This instrument must be installed on a stable table. Do not install this instrument in a location that is unstable or exposed to vibration.</p>
 CAUTION	<p>Federal law restricts the sale of this instrument to sale by or on the order of a physician or practitioner (CFR 801.109(b)(1)).</p>
 CAUTION	<p>When adjusting the Fundus camera head, be careful not to pinch hands with the LCD monitor/touch panel or Base.</p>
 CAUTION	<p>When adjusting the Fundus camera head, be careful not to bring the Ocular lens into contact with the patient's face.</p>

2.2 Symbols and Labels

	Presence of electrical shock hazard. Note: Indicates the risk of electrical shock due to the presence of uninsulated high voltage inside the instrument. Do not remove the instrument cover or parts.
	Circuit breaker
	Type B applied parts. Note: This instrument complies with the specified requirements to provide protection against electrical shock, particularly regarding the allowable patient leakage current.
	Manufacturer Crystalvue Medical Corporation No. 116, Ln. 956, Zhongshan Rd., Taoyuan Dist., Taoyuan City 33072, Taiwan
	Authorized representative in the European Union Medical Device Safety Service (MDSS) GmbH Schiffgraben 4130175 Hannover, Germany
	Serial number
	Catalog number / part number
	European conformity

2.3 Protective Packing Symbols

The protective packing symbols specify the handling requirements and the transportation and storage conditions.

	Fragile, Handle with Care
	Keep Dry
	This end up.
	Relative Humidity
	Temperature
	Do not stack.
	2 Layers Only
	<p>Waste Electrical and Electronic Equipment (WEEE) Recycling Instructions</p> <p>When it has been determined that this instrument is ready for disposal, it must be recycled in accordance with the requirements of the policies and procedures of the corresponding country. Do not dispose of this instrument as general waste.</p>

Product Compliance

	<p>93/42/EEC/M5 Medical Device Directive</p>
	<p>Indicates this equipment contains Type B applied parts</p>
	<p>The RetinaStation is classified as follows:</p> <p>Class I Equipment – Protection against electrical shock.</p> <p>Type B – Degree of protection against electric shock of applied part (chin and forehead rests).</p> <p>Ordinary Equipment (IPX0) – Degree of protection against ingress of liquids (none).</p> <p>Continuous Operation – Mode of operation</p> <p>Electromagnetic Compatibility (EMC): EN 60601-1-2 The RetinaStation device has been tested to comply with the emission and immunity requirements of EN60601-1-2. The RetinaStation is intended for use in an electromagnetic environment where radiated RF disturbances are not beyond the standard defined in EN60601-1-2.</p>

Certification: Under IEC 60601-1

<p>Guidance and manufacturer's declaration – electromagnetic emissions</p>		
<p>The RetinaStation is intended for use in the electromagnetic environment specified below. The customer or the user of the RetinaStation should assure that it is used in such an environment.</p>		
<p>Emissions test</p>	<p>Compliance</p>	<p>Electromagnetic environment – guidance</p>
<p>RF emissions CISPR 11</p>	<p>Group 1</p>	<p>The RetinaStation uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference to nearby electronic equipment.</p>
<p>RF emissions CISPR 11</p>	<p>Class A</p>	<p>The RetinaStation is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.</p>
<p>Harmonic emissions IEC 61000-3-2</p>	<p>Class A</p>	
<p>Voltage fluctuations/flicker emissions IEC 61000-3-3</p>	<p>Complies</p>	

Guidance and manufacturer's declaration – electromagnetic immunity			
The RetinaStation is intended for use in the electromagnetic environment specified below. The customer or the user of the RetinaStation should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Interruptions and voltage variations in power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the RetinaStation requires continued operation during power mains interruptions, it is recommended that the RetinaStation be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the AC mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity

The RetinaStation is intended for use in the electromagnetic environment specified below. The customer or the user of the RetinaStation should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the RetinaStation, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance $d = 1.2 \sqrt{P}$</p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^(a) should be less than the compliance level in each frequency range. ^(b)</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures objects and people.

^(a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the RetinaStation is used exceeds the applicable RF compliance level above, the RetinaStation should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the RetinaStation.

^(b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distances between
portable and mobile RF communications equipment and the RetinaStation**

The RetinaStation is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the RetinaStation can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the RetinaStation as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23


For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

2.4 Product Labels

RetinaStation system label (sample only):

WARNING: Do not connect any unspecified devices to the instrument. Doing so may result in fire or electric shock. For details regarding purchasing accessories, please contact your local Nikon representative or distributor.


2.5 Service Life

The service life of RetinaStation is five years if the specified inspections and maintenance are performed.

3 Description of the Instrument

3.1 Introduction


3.2 Standard Accessories

Item	Description	Quantity
1	User manual	1 pc.
2	AC power cord	1 pc.
3	Dust cover	1 pc.
4	Forehead adaptor	1 pc.
5	Lens cover	1 pc.


4 Introduction to the User Interface

Three tabs are displayed on the screen after the instrument is powered on. These tabs are described in the following sections.


4.1 Managing Patients

Select a patient, add a patient, and modify patient information.


Patient list: Displays all the patients associated with the search result.

Search: This is a function for searching for patients by entering keywords.

- ◆ Search all fields: Enter a keyword such as "World".

- Add patient  : Click this button to add a new patient profile and information.


The 'NEW PATIENT' form contains the following fields:

- *Patient ID (required)
- EMRID
- Name
- Birthday(yyyy-MM-dd)
- Gender (dropdown menu)
- Phone
- Email
- Address
- Comment

Buttons: Close (X), Save (floppy disk icon)

- ◆ Fields marked with a * must be entered to add a new patient to the database.
 - ◆ Click the  button to save the new patient information to the database. The display returns to the main patient information window and the newly added patient is listed in the patient list. Select the newly added patient to display the patient details.
 - ◆ Click the  button to return to the main patient window without saving the patient information.
- Find patient  : Click this button to display the Find patient from DICOM server window.

The 'Find patient from DICOM server' window includes:

- Search conditions: Patient ID, Start Date, End Date, and a SEARCH button.
- Found patient list: A panel showing 'NO PATIENT' with a list icon and a checkmark icon.
- Patient details: A panel showing 'NO PATIENT' with CANCEL and OK buttons.

Labels with arrows point to the 'Found patient list' and 'Patient details' panels.

◆ Search conditions

- Patient ID: Search for patients by patient ID.
- Start date and end date: Search for patients by examination date.

Patient ID

Start Date 2018-10-30

End Date 2018-11-01

SEARCH

◆ Found patient list

- Patient list: This list displays all the found patients.
- Select all and deselect all: Select or deselect all candidate patients.

Candidate patient

No check button means some incorrect data was found.

Not a candidate patient


Click a patient to check the patient's details.

Examination date

Check	Patient Name	Examination Date
<input checked="" type="checkbox"/>	Fundas^Camera^^^	2019-01-11 04:01:56
<input type="checkbox"/>	Pa^Tient^^^	2019-01-11 04:01:06
<input type="checkbox"/>	Pat^lent^^^	2019-01-11 03:01:45


◆ Patient details

- Error message: This is displayed when the found patient contains incorrect data.
- Patient ID: Patient ID of the found patient
- Name: Name of the found patient
- Birthday: Birthday of the found patient
- Gender: Gender of the found patient
- OK button: Submit all candidate patients.
- Cancel button: Close the window.

The screenshot shows a form with the following fields: Patient ID (000003), Name (Pa^Tient^^^), Birthday (2000-10-01), and Gender (Male). At the top, a red error message reads: "Patient ID has been found but patient information is different." Below the form are "CANCEL" and "OK" buttons. A red line points from the error message to the text "Error message" on the right.

All candidate patients are displayed in the patient list after the candidate patients are submitted.


The screenshot shows a patient list interface. At the top is a blue header with a person icon and the word "PATIENT". Below it is a search bar with the text "Search Patient" and "camera". There are two buttons: "SHOW ALL" and "SHOW TODAY". A patient entry is shown with a green refresh icon in a red box, a person icon, the name "Fundas^Camera^ ^^", and "0 image(s)". A red line points from the text "Patient found on DICOM server" to the refresh icon.

NOTE: The patient is created automatically if the patient does not exist on RetinaStation.


■ Patient information

Displaying selected patient information


- ◆ Edit patient information  : Click this button to edit patient information and save comments. Operations are the same as performed for adding a patient.
- ◆ Delete a patient  : Click this button to delete a patient's information and images. When the delete dialog box appears, click [OK] to delete the patient information or [Cancel] to quit.


- ◆ Settings  : Click this button to show the RetinaStation program settings.
 - PACKING: Enable Auto Packing to move the camera head to the packing position when RetinaStation shuts down.


- EXPORT: Enable Auto Export to copy the captured image automatically to the specified folder.

- ARCHIVE: Create a backup file of the RetinaStation database.
 1. Choose the target folder where the backup file will be saved.


2. Click the [ARCHIVE] button to create the backup file.


NOTE: The file system of the target drive should be NTFS.

- **RESTORE:** Restore the RetinaStation database from a backup file.


1. Choose the target folder where the backup file is saved.


2. Choose the backup file from the list and then click the [RESTORE] button to restore the RetinaStation database.


- **DICOM SETTING:** Set the parameters for DICOM communication.

- LANGUAGE: Select the user interface language.


- ABOUT: Display the software version of RetinaStation.

4.2 Capturing Images

Main window for capturing images


Control buttons in the window

- Chin rest: Move the chin rest up and down.
- Z control: Move the camera forward or backward.
- X/Y control: Click the center of the pupil on the screen to perform alignment.
- Default position: Reset the camera to the default position.
- Fixation setting: Select the fixation positions.
- OS/OD selection: Select OD or OS for image capturing.
- Start capture: Click to perform automatic eye alignment and capture an image.
- Mode selection: Select auto alignment mode or manual alignment mode.
- Advanced settings: Display the RetinaStation program settings.
 - ◆ Semi Auto mode: The user switches the eye manually.
 - ◆ Full Auto mode: The RetinaStation switches to the other eye automatically.
 - ◆ Enable manual mode helper: Enable this for help finding a suitable working distance and performing capture processes automatically.


Settings

Mode | Semi Auto | Full Auto
Manual Mode | Enable manual mode helper


- Capture mode control: Display the capture mode setting panel and select either retina (default) or cornea mode.

- Diopter control: Display the diopter setting panel and select the empty lens (default), - lens or + lens.

- Lighting control: Display the lighting setting panel and select the viewing and flash LED levels.


- Dilate selection: Captured images are saved with this setting.

4.2.1 Automatic alignment and focus operation


1. Adjust the chin rest with the   buttons to a suitable position.
2. Click the pupil position on the screen to align the camera.
3. Click the  button to start tracking and capture.

4.2.2 Montage mode


1. Press and hold the Start button for 3 seconds to display the Montage button.


2. Click the Montage button and select the fixation LEDs.

3. Click the [SET] button.
4. Click [OK] to start montage capture.


5. Click the [CONTINUE] button for the next shot.


4.2.3 Manual alignment and manual capture

1. Change the mode selection button to manual mode.
2. Adjust the chin rest with the   buttons to a suitable position.
3. Click the on-screen pupil position or press the     buttons to align the camera.

4. Approach the pupil with the   buttons until the split bar is visible.


5. Align the split bar with the   buttons. (Click on the  button to remove the fixation mask plate. The user can see the live retina video.)

6. Enter the pupil using the   buttons until two spots have appeared.


7. Align the two spots with the       buttons and click the  button to capture an image.


8. Click on the  button to display operation tips.

Manual Mode Tips


4.2.4 Cornea capture mode

1. Add the forehead adaptor to forehead rest.
2. Click the  button and then select the cornea capture mode.


3. Adjust the chin rest with the   buttons to a suitable position.

4. Click the on-screen pupil position or press the     buttons to align the camera.


5. Use the   buttons until the image is clear.

6. Click the  button to capture an image.


4.2.5 Send DICOM images

When the selected patient was found on the DICOM server, the DICOM sending window is displayed before the user leaves the image capturing image window. The user can choose images and send them to the DICOM server.


4.2.6 Review captured images

Display the visit date specified by the capture date. The capture time is also displayed for each image. The operator can review images in this window on a particular date.

■ Multiple selections


1. Press and hold on an image thumbnail to enter the multiple selection mode.
2. Click an image thumbnail to select that image for batch deletion or exporting.

DICOM Montage Printer Compare Export Delete Cancel selection


■ Delete an image

Click the  button. A warning dialog appears. Click [OK] to delete the selected image. Be careful, as a deleted image cannot be restored!


■ Montage

1. Click the  button to display the montage dialog box.

2. Click the  button to display a visit date list for the selected patient.

Visit Date list

Click on an image to assign it to the candidate list.

Candidate list


Start montage

3. Click the  button to display the montage result.

- Export an image

Click the  button to display a warning dialog box. Then choose the destination folder and image compression format. Check Export ID to include the patient ID in the file name.


- Export all images of a patient


Press and hold the patient name in the patient list. The export button appears. Click the  button and finish export options,

■ Compare images

- ◆ Click the  button to display the comparison dialog box.


- ◆ Date list switch: Click the  button to display a visit date list for the selected patient.

Visit date list

Click on an image to assign it to the left or right frame.


Left frame

Right frame


- ◆ Reset: Click the  button to reset all image settings to the default settings.
- ◆ Redfree: Click the  button to show the redfree image.

■ Printer

- ◆ Click the  button to display the printer dialog box. Click [PRINT] to print the selected image.

- ◆ To print two images, select the images and click the  button, then click [PRINT] in the printer dialog box.


■ DICOM

- Click the  button to send the selected image(s) to the DICOM server.


4.2.7 Photo viewer

Click an image in the image list to display the photo viewer. The photo viewer has the following functions.


- Image information
ID: Patient ID
Eye: Captured eye
Mode: Capture mode
Time: Date and time of capture

- Close: Exit the photo viewer.
- Contrast: Adjust the selected image. The  button increases the contrast level and the  button decreases the contrast level.


- Brightness: Adjust the selected image. The  button increases the brightness level and the  button decreases the brightness level.

- Negative: Do negative processing for the selected image.
- Redfree: Remove the red channel and convert the selected image to grayscale.

- RGB separation: Do RGB channels separation and display in color or grayscale mode. The  button is grayscale mode and the  button is color mode.

- CD ratio: Perform disc and cup measurement. The following describes how to measure the cup to disc ratio.
 1. Use the pinch gesture to zoom in/out on the region of interest.
 2. Click the CD ratio button and click the  button.
 3. Tap the border of the disc until the blue line fits it.

4. Click the CD ratio button and click the  button.
5. Tap the border of the cup until the yellow line fits it.

6. Click the CD ratio button and click the  button to save the measurement result. The measurement is displayed on the top right side of the photo viewer.

Click the  button to clear the measurement result.

	ID	000001
	Eye	OD
	Mode	Retina
	Time	2018-11-07 13:50:46
Cup to disc horizontal ratio	CD H	0.7
Cup to disc vertical ratio	CD V	0.6

- Export: Export the current image.
- Reset: Reset the measurement, scaling ratio, image position and processing level etc.
- Back/Next: Choose the next or previous image.

5 Specifications

5.1 Product Specifications

Function	Value/Type	Remark
Fundus image		Non-mydratiac, color image
Field of view	>= 45 degrees	
Illumination for retina image (capture)	White LED	Flashes for capturing images
Cornea image (capture)	White LED	Flashes for capturing images
Illumination during alignment to patient's retina	NIR LED	Central wavelength in the range of 735 – 850 nm
Focus diopter adjustment range	-15 D to +10 D	Without compensation lens
	-30 D to -10 D or +5 D to + 30 D	With compensation lens
Minimum pupil size	4 mm	
Focus adjustment	Auto / manual	Split-image technique
Image sensor	CMOS 12 megapixel	
Z-ranging (working distance)	2 fiber dots	
Working distance	25 mm from lens to cornea	Accuracy: +/- 0.5 mm
Fixation	Internal	10 points

General

Function	Value/Type	Remark
Alignment	Fully automatic 3D tracking	
Alignment mode	Full auto / auto/ manual	
Chin rest	Motorized	
Interface	USB 2.0 port, LAN, HDMI	
Input/output format	Image format: JPEG, PNG, BMP, DICOM	
Display	10.1" LCD monitor, touch panel	
Operation range	Front / back: 40 mm	
	Left / right: 90 mm	
	Up / down: 30 mm	
Chin rest range	Up / down: 70 mm	
Dimensions	Depth: 495 mm	
	Width: 288 mm	
	Height: 495 mm	
Weight	17 kg	

Note 1: The USB interface is used for connecting a USB mass storage device.

5.2 Environmental Conditions

5.2.1 Operating Conditions:

- Temperature: 10°C to 35°C
- Humidity: 30% to 90% RH
- Atmospheric pressure: 800 – 1060 hPa

5.2.2 Storage Conditions:

- Temperature: -10 – 55°C
- Relative humidity: 10 – 95% RH
- Atmospheric pressure: 700 – 1060 hPa

5.2.3 Transport Conditions:

- Temperature: -40 – 70°C
- Relative humidity: 10 – 95% RH
- Vibration, sinusoidal: 10 – 500 Hz, 0.5 G
- Shock: 1/2 Sine wave, 6 msec, 30 G peak (packaged)
- Bump: 1/2 Sine wave, 6 msec, 10 G peak (packaged)

5.3 Electric Rating

- Source voltage: AC 100 – 240 V
- Frequency: 50 – 60 Hz
- Power input: < 150 VA

6 Maintenance

6.1 Cleaning the Lens

It is recommended that the ocular lens of the RetinaStation is cleaned on weekly basis or regularly as needed.

- Items required for cleaning the ocular lens:
 - a) Diluted acetone or lens cleaning solution
 - b) Lens cleaning paper
- Method:

Wet a sheet of lens paper with some cleaning solution and wipe the ocular lens with one pass in one direction. Discard the used sheet of lens paper. Use a new sheet for each repeat cleaning until the ocular lens is clean.

6.2 Chin Rest and Forehead Rest

Using a cleaning cloth dampened with a disinfectant solution or a wet isopropyl alcohol cleaning paper pad, wipe the chin rest and forehead rest before and after each use.

If chin rest paper is used, remove one piece for each patient.

6.3 LCD Monitor/Touch Panel

First, turn off the power and, using a soft cleaning cloth, lightly wipe the exterior of the LCD display. Do not press the LCD monitor with a hard tipped object. Doing so may scratch or damage the LCD monitor, resulting in malfunction.

6.4 Changing the Fuse

1. Check that the power is turned off and the power cable is unplugged.
2. Tilt the body slowly so that the AC Power Entry on the base can be seen.
3. Remove the fuse holder from the AC Power Entry.
4. Replace the fuse with the correct part.

Fuse specification:


- Rating: 1.6 A / 250 VAC
 - Package: 5 mm x 20 mm
 - Type: Slow blow
5. Insert the fuse holder back into the AC Power Entry.

6.5 Cleaning the enclosure


If the enclosure needs cleaning, wipe the soiled area using a cleaning cloth dampened with a disinfectant solution or a cleaning paper pad dampened with isopropyl alcohol.

7 Updating the Software


1. Plug the USB drive into a USB port and copy the RetinaStation installation package file to the Windows desktop.


2. Double click the RetinaStation installation package file. The following installation dialog box appears.

3. Click the [Install] button and wait for installation to complete.

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